

Report info**Participants**

61 participants from 15 countries.

Report info

Qualitative and quantitative results are presented in the same report. Always compare your results to those obtained with the same test kit.

The qualitative results are divided into groups according to the method stated by the laboratory and presented in laboratory-specific tables. Accepted results are marked with green colour and laboratory's own result with a black radio button .

The quantitative results from the positive specimens are presented as histograms. Quantitative results reported with < or > -signs are not processed. Results reported as ratios, titres, absorbances and without a unit are manually removed and not statistically processed.

Assigned values (target values) are mean of the results where results deviating more than +/- 3*standard deviation from the median are removed. The standard uncertainty (u) of the assigned value is reported as standard error of the mean (SEM). Additionally, if the measurement uncertainty of the target value is large an automatic text is printed on the report: "The uncertainty of the assigned value is not negligible, and evaluations could be affected". In case the client's result is the only one in the method group, no assigned value will be calculated, no target area shown, and no statistics calculated.

In participant specific report your own result is shown with an orange dot. The target area is presented as a yellow area in the picture. In the history graphs you are able to see your performance graphically with both against the assigned value (xpt) and the z-score area of -2 ... +2.

If you have not reported any results you will get a note: "You have not responded in time, only global report is available".

For information on report interpretation and performance evaluation, please see the "EQAS Interpretation guidelines" in LabScala User instructions. In case you have any questions regarding the reports, please contact the EOA Coordinator.